LS810-B

TENSIÓMETRO DE MUÑECA INALÁMBRICO
WIRELESS WRIST BLOOD PRESSURE MONITOR
TENSIÓMETRO DE PULSO SEM FIO
TENSIO MÈTRE POUR POIGNET SANS FILS

- ES GUÍA DE USUARIO.....2
- EN USER GUIDE.....13
- PT MANUAL DE INSTRUÇÕES.....24
- FR GUIDE DE L'UTILISATEUR.....35

LA CAJA CONTIENE:

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THE BOX CONTAINS

An instruction manual, IVBPW wrist blood pressure monitor and USB/micro USB charging cable

SOME TECHNICAL DETAILS

Bluetooth Smart 4.0 (low energy)

Accuracy Pressure: 5°C-40°C within ± 3 mmHg (0.4kPa)

Pulse value: $\pm 5\%$

Blood pressure measurement range: 40-230 mmHg

Pulse rate measurement range: 40 - 199 beats/minute

Temperature/ humidity for use: +5° C a +40° C / max: $\leq 85\%$ RH

Temperature/humidity for storage: -20° C a +60° C / max: 10%RH-93%RH

On-line and off-line operation

INTENDED USE

The IVY HEALTH IVBPW blood pressure device allows you or your doctor to take your blood pressure using non-invasive means. All you have to do to charge the device's internal battery is connect the cable included to a charging socket or PC USB port. It's that simple: the blood pressure device is then ready to take your blood pressure and pulse rate.

SAFETY WARNINGS

Do not use this device for self-medication on the basis of its measurements or for any therapeutic measure without consulting your doctor.

Contraindications:

1. This blood pressure monitor should not be used if you:
 - a. Suffer from severe arrhythmia, pre-eclampsia, ventricular and/or arterial occlusive disease.
 - b. Are undergoing intravascular therapy
 - c. Have venous access
 - d. Have had a mastectomy
 - e. Are pregnant
2. Do not allow its use by children
3. If you are suffering from any illness, consult your doctor before using this blood pressure monitor.



Safety warnings on the battery and other parts of the device:

1. Never open the blood pressure device or dismantle the batteries
2. Keep the batteries out of reach of children. Do not throw them in fire or short circuit them as they could explode.
3. Do not throw the batteries or any of the device components in ordinary rubbish. Take them to specialised collection centres or put them in a proper container if there is one in your area.
4. If someone accidentally swallows a battery, please try to find medical assistance immediately.

5. If the battery level is inadequate, it can give rise to incorrect measurements or communication errors with the IvyHealth application.
6. Do not use the blood pressure device if you are allergic to any of its components.
7. Do not use the blood pressure device for purposes other than those indicated in this manual.

Cleaning the device:

1. Use a damp cloth with a soap solution, but make sure no liquid gets into the device. The blood pressure device should not be used if it is wet.
2. Use a damp cloth with a soap solution, but make sure no liquid gets into the device. The blood pressure device should not be used if it is wet.
3. You are recommended to disinfect the inner part of the cuff in contact with the skin with a soft cloth moistened with ethyl alcohol every once in a while.
4. Do not inflate it until it has been correctly positioned on your arm.
5. Do not expose it directly to the sun, extreme humidity and/or dust.


WHAT CAN I DO WITH THIS BLOOD PRESSURE DEVICE?

Can it work without the application and my mobile phone?

Yes, you can use the blood pressure device without the free IvyHealth app and without connection to a mobile phone or tablet. You can see the results on the device's screen. However, if you are connected to a mobile phone with the IvyHealth app open, you will have a more complete view of your blood pressure and other values.

You are recommended to take the readings at the same time and under the same conditions in order to be able to compare results. And it is always better to take your blood pressure before meals and after resting for a few minutes.

The blood pressure device has an internal memory where all the measurements taken, together with the date on which they are made, are automatically saved.

To adjust the blood pressure device's default date, keep the SET key pressed with the monitor switched off until the device switches on. The time displays in HH:MM format with the time field flashing, and the date in MM/DD/YY format. Press the MEM key as often as necessary to obtain the correct time value. Once you have reached it, press SET and the next field to be modified will start flashing. Continue with the process until you have gone through all the time and date fields. Once you have finished, the word DONE displays on the screen implying that the time and date have been modified.

To see previous measurements on the screen, press MEM with the blood pressure device switched off. This will switch on the blood pressure monitor and display the measurement with the most recent date. The screen will display the values and the time and date on which they were taken. Press SET to access the previous record and MEM the subsequent one.

To delete all the records and reset the blood pressure monitor to factory settings, switch off the device and press MEM and START/STOP at the same time. Sh1 P displays on the screen. Then keep the SET key pressed until the blood pressure device switches off. To switch it on again for the first time following memory clear, keep the SET and START/STOP keys pressed. Remember that you will need to reset the time and date.

What if I download the application?

Download the free IvyHealth application from the Apple Store or Google Play store. When you first open the application, you must register and set up your user name. It is very easy if you follow the instructions on your mobile phone's screen.

The IvyHealth application can be automatically paired with your recently acquired blood pressure device.

1. When you first open the application you'll find the registration screen. Once you are registered you will enter the main menu, where you'll be able to choose the device you'll be working with.

2. To make a measurement, press the icon  on the main screen and "New measurement". You will access the measurement screen and can now position the blood pressure monitor on your arm and follow the instructions on the screen.
3. Press the START/STOP key on the blood pressure device.
4. The blood pressure monitor screen will display the provisional data on systolic pressure while the device is electronically inflated. Your telephone or tablet will automatically detect the measurement and the collected data will be shown in the screen.

5. Wait until the measurement is completed.
6. Once the transmission has been completed, you can see the following data (*) on your mobile phone screen:
 - a. Your systolic pressure
 - b. Your diastolic pressure
 - c. Your pulse rate
 - d. IHB (Irregular Heart Beat), an indicator of an irregular heart beat or arrhythmia.
 - e. Once the measurement is finished, click on the "Report" button and you will see a highly intuitive report that tells you if the readings made are at the required level or if your blood pressure is too high or too low, in accordance with the readings laid down by the WHO.
7. Press to save and your first reading is recorded for the future. You can consult your phone whenever you wish. You will be protected by your own password.
8. Subsequent measurements will be kept in your mobile phone whenever you measure your blood pressure with the app open and you follow the previous steps.

What else can I do with the application?

If you use the blood pressure monitor together with the free IvyHealth application, you can obtain the readings for all the body values listed in the previous paragraph. However, the application also allows you to:

1. Consult your statistics. Just click on the icon  and then access to the "History" section.
2. Save and display the results of as many users as you wish. Every time somebody makes a measurement from their mobile phone, and if they are registered as a user, they can record their results and protect them with a password.

THE MOST COMMON PROBLEMS AND THEIR SOLUTION

The blood pressure monitor does not respond to any button: Make sure the battery is not flat. If this is the case, charge the blood pressure monitor. If the problem persists, keep the SET and START/STOP buttons pressed until the blood pressure device switches on.

The blood pressure monitor displays the message:

 + Lo / ERROR The battery is wearing out. Please, charge the monitor.

ERROR 1 / ERROR 2 / ERROR 3 / ERROR 21 This is due to the fact that insufficient pressure has been detected in the inflation process or the sleeve is not adjusted properly. Try removing and readjusting the monitor, as shown in the drawing of the cuff.

ERROR 5 / ERROR 6 System error. Please, try again. If the problem persists, contact our customer service department for further assistance. Refer to the warranty for contact information and return instructions.

ERROR 10 It is due to a problem in the inflation process. This could be caused by some movement or adjustment once measurement has started.

ERROR 20 The pulse rate is not detected. Make sure there isn't any cloth between your wrist and the monitor and tray again.

In addition to this

The IVY HEALTH package of health complements is made up of different devices affording you greater control over your health.

Check out our latest products
for health care in this link.
www.ivyhealth.es

LEGAL NOTICE

Declaration of Conformity

FONEXION SPAIN S.A. declares that the IVPBW digital wireless blood pressure device is in compliance with European standards EN 1060-1 and EN 1060-4 as well as with the essential requirements and other relevant provisions of the European Directive MDD93/42/CE, as well as with Directive R&TTE 1999/5/CE. You can download a copy of the complete declaration at: <http://www.fonexion.com/doc>.

Guarantee

The IVPBW digital wireless blood pressure device has a two-year guarantee from the date on the sales receipt. The sales receipt is your guarantee, so keep it. Nevertheless, this guarantee does not cover any fault caused by accidents, improper use, wear, negligence, lightning, modifications made to the equipment, or any attempt at adjustment or repair in the home or by unauthorised agents.

Bluetooth

The Bluetooth brand and its logos are registered by Bluetooth SIG, Inc. and use of said brands by IVY HEALTH is made under licence. Other brands and trade names belong to their respective owners.

STANDARDS AND DOCUMENTS

Safety Information:

	Symbol for "THE OPERATION GUIDE MUST BE READ"		Symbol for "TYPE BF APPLIED PARTS"
CE 123	Symbol for "COMPLIES WITH MDD 93/42/EEC REQUIREMENTS"		Symbol for "ENVIRONMENT PROTECTION - Electrical waste products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice"
	Symbol for "MANUFACTURER"		
SN	Symbol for "SERIAL NUMBER"		Symbol for "DIRECT CURRENT"
	The Bluetooth Combination Mark		Symbol for "Authorised Representative in the European Community"
	Symbol for "MANUFACTURE DATE"		Symbol for "Including RF transmitter"

Specifications:

Power supply	3.7V 420mAh Built-in rechargeable lithium-ion battery, 5V \equiv 1A USB cable
Display mode	Digital LCD V.A.46.5mmx36.5mm
Measurement mode	Oscillographic testing mode
Measurement range	Rated cuff pressure: 0mmHg-300mmHg(0kPa-40kPa) Measurement pressure: 40mmHg-230mmHg (5.3kPa-30.7kPa) pulse value:(40-199)beat/minute
Accuracy	Pressure: 5°C-40°C within \pm 3mmHg(0.4kPa) Pulse value: \pm 5%
Normal working condition	Temperature:5°C to 40°C Relative humidity: \leq 85%RH Atmospheric pressure: 86kPa to 106kPa
Storage & transportation condition	Temperature:-20°C to 60°C Relative humidity: 10%RH to 93%RH Atmospheric pressure: 50kPa to 106kPa
Measurement perimeter of the wrist	About 13.5cm-21.5cm

Weight	Approx.110g
External dimensions	Approx.79.8mm×72.5mm×13.2mm
Attachment	USB cable, user manual
Mode of operation	Continuous operation
Degree of protection	Type BF applied part
Protection against ingress of water	IP22, It means the device could protect against solid strange objects of 12.5 mm and greater, and against vertically falling water drops when ENCLOSURE tilted up to 15°
Software version	V01
Device classification	Battery Powered Mode: Internally Powered ME Equipment AC Adaptor charged Mode: Class II ME Equipment

Standards:

Title	Code
Medical Device Quality System for the regulation requirement	EN ISO 13485:2012
Medical device risk management to medical devices application	EN ISO14971:2012
Symbol for the label of medical devices	EN 980:2008
Term, symbol and information of medical device— information of medical device manufacturer offering	EN 1041:2008
Medical electrical equipment – Part 1: General requirements for safety	EN 60601-1: 2005/AC:2012
Medical electrical equipment – Part 1-2: General requirements for safety – Collateral standard: Electromagnetic compatibility – Requirements and tests	EN 60601-1-2:2007
Clinical investigation of medical devices for human subjects. General requirements	EN ISO 14155-1:2009
Clinical investigation of medical devices for human subjects. Clinical investigation plans	EN ISO 14155-2:2009
Medical device software – Software life-cycle processes	EN 62304:2006/AC: 2008
Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral standard: Usability	EN 60601-1-6:2010
Medical electrical equipment. Part 1-11:General requirements for basic safety and essential performance. Collateral standard:Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment.	EN 60601-1-11:2010

Medical devices – Application of usability engineering to medical devices	EN 62366: 2008
Council Directive 93/42/EEC of 14 June 1993 concerning medical devices OJ L 169 of 12 July 1993 (including Directive 2007/47/EC)	MDD 93/42/EEC
Non-invasive sphygmomanometers —Part 1: General requirements	EN 1060-1
Non-invasive sphygmomanometers.Part 3:Supplementary requirements for electro-mechanical blood pressure measuring systems	EN 1060-3
Secondary cells and batteries containing alkaline or other non-acid electrolytes – Safety requirements for portable sealed secondary cells, and for batteries made from them, for use in portable applications	IEC 62133:2012
User Guide	Versión 1.0

EMC GUIDANCE:

**Guidance and MANUFACTURER'S declaration – ELECTROMAGNETIC EMISSIONS-
for all ME EQUIPMENT and ME SYSTEMS:**

Guidance and manufacturer's declaration – electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 2	The device must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Guidance and MANUFACTURER'S declaration – electromagnetic IMMUNITY – for all ME EQUIPMENT and ME SYSTEMS:

Guidance and manufacturer's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 s	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE UT is the a.c. mains voltage prior to application of the test level.			

Guidance and MANUFACTURER's declaration – electromagnetic IMMUNITY –
for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING:

Guidance and manufacturer's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 Vrms</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance: $d = 1.167 \sqrt{P}$ $d = 1.167 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.333 \sqrt{P}$ 800 MHz to 2,5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a) should be less than the compliance level in each frequency range b).</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>
<p>NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>a. Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.</p> <p>b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.</p>			

Recommended separation distances between portable and mobile RF communications equipment and the ME EQUIPMENT or ME SYSTEM –for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING:

Recommended separation distances between portable and mobile RF communications equipment and the device.			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.167 \sqrt{P}$	80 MHz to 800 MHz $d = 1.167 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.333 \sqrt{P}$
0.01	0.117	0.117	0.233
0.1	0.369	0.369	0.738
1	0.369	0.369	2.333
10	3.690	3.690	7.378
100	11.67	11.67	23.33
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80MHz and 800MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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FONEXION SPAIN S.A. GRAN VÍA 45, 1º PLANTA DPT.114. 48011 BILBAO
Fabricado para FONEXION SPAIN S.A. VAT:ESA48937239 Gran Vía 45, 1-114, 48001 Bilbao (SPAIN) · www.fonexion.com
MDSS · Medical Device Safety Service GmbH · Schiffgraben 41 · 30175 Hannover · Germany
FONEXION SPAIN S.A. Gran Vía 45, 1-114, 48001 Bilbao (SPAIN) · www.fonexion.com